



ProFoldin

10 Technology Drive, Suite 40, Number 188
Hudson, MA 01749-2791 USA
Tel: (508) 735-2539 FAX: (508) 845-9258
www.profoldin.com
info@profoldin.com


INSTRUCTIONS

ProFoldin Easy Protein Assay Reagent

CATALOG NUMBER EPA001

INTRODUCTION

The Easy Protein Assay Reagent provides a quick and simple method to measure protein concentrations. The assay is based on protein - Coomassie Blue binding and its light absorbance 595 nm. The assay is quick and easy: simply mix the diluted reagent with 1/10 volume of the protein solution and read the light absorbance at 595 nm. Standard clear 96-well plates or cuvettes can be used for absorbance reading.



The Easy Protein Assay Reagent (Catalog number EPA001) includes 100 ml of 5 x reagent. It is for 2500 assays using 96-well plate or 500 assays using 1 ml-cuvette.

PROTOCOL

- (1) Diluted 1 ml of the reagent with 4 ml of water.
 - (2) Mix 200 µl of the diluted reagent with 20 µl of the protein solution for measurement using 96-well plates. For measurements using 1 ml- cuvettes, mix 1 ml of the diluted reagent with 0.1 ml of the protein solution.
 - (3) Read the light absorbance at 595 nm (A_{595}).
-



ProFoldin

10 Technology Drive, Suite 40, Number 188

Hudson, MA 01749-2791 USA

Tel: (508) 735-2539 FAX: (508) 845-9258

www.profoldin.com

info@profoldin.com

INSTRUCTIONS

Data Analysis

Plot the A_{595} and the protein concentration [**Protein**] to generate the linear standard curve.

$$A_{595} = a [\text{Protein}] + b$$

Where the A_{595} values are from experimental data, the **a** and **b** values are from the linear fitting between the A_{595} values and the protein concentrations.

UNKNOWN SAMPLES

Follow the same procedure to measure the A_{595} values from the unknown samples. Calculate the protein concentrations in the unknown samples using the A_{595} values from the unknown samples and the **a** and **b** values from the standard curve.

$$[\text{Protein}] = (A_{595} - b) / a$$